



2022

In partnership with



Test kits

Detection of allergens, food pathogens, Mycotoxins, GMOs,
Antibiotic residues in milk and other Neogen products



1-888-8BIOVET

(824-6838)

order@biovet-inc.com

Summary

Summary	2
Term and conditions	3
About Biovet	3
Allergens	4
Reveal.....	4
Alert.....	5
Veratox.....	6
Mycotoxins.....	7
Reveal.....	7
Veratox.....	8
Others*	8
GMOs and ruminants.....	9
Water	9
Colitag	9
Sanitation	10
Accuclean Advanced	10
AccuPoint Advanced Next Generation.....	11
Others	11
Instruments & material.....	12
Readers	12
Pipettes	12
Others	13
Other products available	15
Pathogens (ANSR)	15
Pathogens (Reveal)	15

MKT-036-EN-V.7 — Revised on 2022-08-02

Term and conditions

A freight charge of \$16.00 is applicable for orders for equipment under \$500.00.

Orders received before Monday 10 am AM will be delivered in the same week. Orders received past this deadline will be deliverable the following week.

Send your order to: order@biovet-inc.com

Effective Date June 20, 2022.

About Biovet

Biovet is proud to offer you unique, quality and innovative products, in collaboration with its partner Neogen. These are powerful tools allowing you to do your tests in-house, for faster results.

Since 1991, Biovet has been operating state-of-the-art laboratories and carrying out microbiological analyses, detection of pathogens, allergens, mycotoxins and GMOs, as well as in the fields of molecular biology and serology (ELISA).

Our laboratory is accredited by the Standards Council of Canada (SCC) for microbiological testing of foods and edibles (see our accreditation scope on our website).

We are multiplying the actions in order to get closer to you with now 2 DIAGNOSTIC SERVICES LABORATORIES in Quebec: SAINT-HYACINTHE and QUEBEC CITY.

Also, our customers appreciate our HIGH AVAILABILITY, since our laboratories are open 6 days and operate 7 days a week.

In addition, Biovet has an efficient collection service that ensures the transport of samples in Quebec as soon as possible, even in the regions.

For further information, feel free to contact your Sales representative.



Allergens

Code	Description
TRD-1017	Allergen Extraction Kit
TRD-1010	Allergen Environmental Swabbing Kit (100 swabs)

Reveal

Neogen's Reveal® and Reveal 3D food allergen detection kits are simple, quick and easy for detecting multiple allergens on work surfaces, ingredients and finished products.



Code	Description	Number of tests
	Reveal 3D	
TRD-1170	Reveal 3D Peanut *	10
TRD-1194	Reveal 3D Almond *	10
TRD-1211	Reveal 3D Crustacea	10
TRD-1186	Reveal 3D Gliadin	20
TRD-1411	Reveal 3D Gluten (New) *†	10
TRD-1221	Reveal 3D Total Milk	10
TRD-1228	Reveal 3D Mustard	10
TRD-1193	Reveal 3D Hazelnut *	10
TRD-1247	Reveal 3D Multi Treenut	10
TRD-1185	Reveal 3D Egg *	10
TRD-1187	Reveal 3D Soy *	10
TRD-1284	Reveal 3D Sesame	10
TRD-1175	Reveal 3D Coconut *	10
	Equipment needed but sold separately from the kit	
TRD-1410	* Reveal 3D Food Buffer	
TRD-1409	† Reveal 3D Gluten Food Buffer	
TRD-1408	* Tubes 50 ml graduated polypropylen (25)	

* Available to analyze on food and ingredients

Code	Description	Number of tests
	Reveal	
TRD-1058	Reveal Histamine	24

Alert

Alert® food allergen detection kits are qualitative ELISA-like screening tests that compare up to 5 samples with a control sample.



Code	Description
TRD-1095	Alert Histamine
TRD-1018	Alert Total Milk
TRD-1013	Alert Egg
TRD-1116	Alert Sulfite (shrimp only)

Veratox

Veratox® food allergen detection kits are quantitative ELISA-like tests that can simultaneously measure up to 19 samples based on a curve consisting of five controls.



Code	Description
TRD-1166	Veratox Allergen Starter Kit
TRD-1016	Veratox VIP for Peanut *
TRD-1014	Veratox Peanut
TRD-1097	Veratox Almond
TRD-1215	Veratox Crustacea
TRD-1196	Veratox Gliadin R5
TRD-1197	Gliadin cocktail solution 25ml
TRD-1205	Gliadin cocktail solution 100ml
TRD-1198	Gliadin cocktail solution 500ml
TRD-1027	Veratox Histamine
TRD-1009	Veratox Total Milk
TRD-1252	Veratox Mustard
TRD-1098	Veratox Hazelnut
TRD-1008	Veratox Egg
TRD-1184	Veratox Soy
TRD-1183	Veratox Sesame

* = New

Mycotoxins

Reveal

The Reveal® and Reveal Q+ test kits are quantitative tests for the analysis of mycotoxins in grains.



Code	Description	Number of tests
TRD-1300	Reveal Q MAX DON	25
TRD-1281	Reveal Q+ DON	25
TRD-1263	Reveal Q MAX Ochratoxin	25
TRD-1299	Reveal Q+ Ochratoxin	25
TRD-1413	Reveal Q MAX T2/HT-2	25
TRD-1288	Reveal Q+ T2-HT2	25
TRD-1266	Reveal Q MAX Zearalenone	25
TRD-1264	Reveal Q+ Zearalenone	25
TRD-1265	Reveal Q MAX Aflatoxine	25
TRD-1207	Reveal Q+ Aflatoxine	25

Mycotoxins

Veratox

Veratox® mycotoxin test kits are quantitative ELISA-like tests that can be assayed simultaneously up to 19 samples based on a curve consisting of five controls.



Code	Description
TRD-1037	Veratox Ochratoxine
TRD-1260	Veratox Ochratoxine grain
TRD-1029	Veratox T-2/HT2
TRD-1002	Veratox Zearalenone
TRD-1004	Veratox Aflatoxine HS
TRD-1028	Veratox Aflatoxine
TRD-1001	Veratox DON 5/5
TRD-1092	Veratox DON HS
TRD-1109	Veratox DON 2/3
TRD-1038	Veratox Fumonisine
TRD-1128	Veratox Fumonisine 5/10

Others*

Code	Description
TRD-1087	Rock-it for mycotoxines
TRD-1287	DON Check Sample (MRM), Low, 100g
TRD-1119	DON Check Sample (MRM), Medium, 100 g
TRD-1110	NeoColumn for Aflatoxine (large) - 50 c

* Check Sample are available for other toxins, for more information, contact your Sales representative

GMOs and ruminants

Code	Description	Number of tests
TRD-1012	Reveal for CP4 (Roundup Ready)	50
TRD-1034	Reveal for Ruminant in Feed	25
TRD-1195	Reveal for Ruminant in MBM	25

Water

Colitag

Code	Description
TRD-1089	Colitag water tes kit P/A – (20)
TRD-1080	Colitag water tes kit P/A – (100)
TRD-1213	Colitag Comparator



Sanitation

Accuclean Advanced

AccuClean® Advanced is a simple test to detect protein residues from liquid or solid food products that have been in contact with the surface to be analyzed. This test is completely autonomous, no instruments are required for reading and it requires only minimal training.



Code	Description
TRD-1171	Accuclean Advanced – Detects ATP (100)

AccuPoint Advanced Next Generation

NEOGEN® New AccuPoint® Advanced Next Generation (NG) is a complete sanitation monitoring system comprising a handheld instrument, innovative samplers, and data manager software, designed to verify your cleaning processes by measuring adenosine triphosphate (ATP). After cleaning your environment, measuring the amount of ATP present on a surface or in liquids allows you to verify your cleaning processes have been effective.



Code	Description
TRD-1311	Accupoint Advance Next Generation
TRD-1242	AccuPoint Advanced Sampler Surface - (100)
TRD-1243	AccuPoint Advanced Sampler Water - (100)
TRD-1244	AccuPoint Advanced Sampler Access - (100)

Others

To check the availability and price of these products, contact your Sales representative.

Code	Description
9950	AccuPoint Advanced ALP (50 Samplers+ 5 controls)
9971	AccuPoint Advanced ALP (5 controls)

ALP : alkaline phosphatase

Instruments & material



Readers

Code	Description
TRD-1042	Neogen Stat-Fax 4700 Microwell Reader
TRD-1292	Raptor Reader (Reader + software)

Pipettes

Code	Description
TRD-1033	Pipettor 100 μ L fixed
TRD-1043	Pipettor 50-300 μ L, 12 Channels
TRD-1067	Pipettor 100 μ L fixed
TRD-1235	Pipettor 100 μ L fixed for Q+ (no overflow)
TRD-335	Pipette Tips (1000)
TRD-1216	Pipettor 500 μ L fixed
TRD-1234	Pipettor 1000 μ L fixed
TRD-1402	Pipettor 400 μ L fixed
TRD-1088	Pipettor 20-200 μ L
TRD-1214	Pipette Tips 100 μ L - 1ml (1000)
TRD-1414	Pipettor 250 μ L fixed
TRD-1401	Pipettes 400 μ L exact volume for Raptor (50)
TRD-1236	Sterile Pipettor 10ml (100)

Others

Code	Description
TRD-1031	Filter Syringes (100)
TRD-1032	Well Holder
TRD-1045	Plastic Wash Bottle, 500 ml (package of 4)
TRD-1047	Incubator with Thermometer
TRD-1049	Microcentrifuge Tube 2 ml (500)
TRD-1056	Whirl-Pak Filter Bags 24oz (250)
TRD-1059	Shaker Water Bath
TRD-1066	Agri-Grind
TRD-1068	Whatman Filter # 1 (100)
TRD-1069	Thermometer
TRD-1070	Graduated Plastic Cylinder, 250 ml
TRD-1071	Plastic Wash Bottle, 250 ml
TRD-1072	Nalgene Bottle 1 litre
TRD-1073	Whatman Filter # 4 11cm
TRD-1091	Electronic Digital Scale
TRD-1094	Premoistened Sampling Sponges
TRD-1096	Timer 3-time
TRD-1101	Allergen Wash Buffer Concentrate
TRD-1105	Tubes, 13 x 100 mm (1000)
TRD-1111	Tubes, 12 x 75 mm (1000)
TRD-1135	Graduated Plastic Cylinder, 50 ml
TRD-1136	Membrane Filters, ISO-GRID, 1600-Grid (box of 100)
TRD-1140	ISO-GRID Filtration Unit
TRD-1142	Blender Jar for mixing slides and lid
TRD-1145	Filtration 6-places
TRD-1147	Vacuum Tubing Assembly
TRD-1148	ISO-GRID Filtration Unit Clamp
TRD-1152	Blender, Oster
TRD-1154	Reagent Boats (3/pack)
TRD-1164	Tubes w/caps 3.5ml (100)
TRD-1165	Amber glass bottle 500ml
TRD-1190	Filtration ramp
TRD-1191	Environmental Sampling Kit
TRD-1202	Pipette Pump 10ml
TRD-1212	Graduated glass Cylinder 50 ml
TRD-1238	Buffer bottle (72/cs)
TRD-1289	Thermal 2 1/4 Paper

Instruments & material / media

Code	Description
TRD-1290	Rack for Reveal Sample Cups
TRD-1400	Raptor Vortex Cartridges (50)
TRD-1025	Fax Printer Paper
TRD-1039	Fax Printer Paper Stat-Fax (10/pk)
TRD-316	Jars and lids for Raptor (100)
TRD-1081	Sample Containers with sodium thiosulfate 100'

Other products available

To check the availability and price of these products, contact your Sales representative directly.

Pathogens (ANSR)

Code	Description
9837	ANSR Complete System
9822	ANSR for E. coli - 96 tests
9815	ANSR for E. coli Enrichment Broth, 500g
9815B	ANSR for E. coli Enrichment Broth, 2 kg
9873	ANSR for Listeria Right No- 96 tests
9824	ANSR for Listeria mono - 96 tests
9814	ANSR Listeria Enrichment Broth 1, 500g (environmental sample)
9880	ANSR Listeria LESS PLUS media, 500g (products and environment)
9870	ANSR for Salmonella - 96 tests
	For Enrichment Broth (Salmonella), see your Sales representative

Pathogens (Reveal)

Neogen's Reveal® and Reveal 2.0 test kits contain everything you need to perform the analysis; including the enrichment medium to detect different pathogens.



Code	Description
9750	Reveal E. coli 0157:H7, 8 h, SC - 20 tests
9751	Reveal E. coli 0157:H7, 20 h, SC - 20 test
9714	Reveal 2.0 E. coli 0157:H7 SC – 20 tests (sample of 65g)
9803	Reveal 2.0 for Salmonella , SC avec RV - 20 tests
9805	Reveal 2.2 for Salmonella , SC avec RV et M-broth - 20 tests
	Reveal 2.0, Listeria One step CS for food - 20 tests
	Reveal 2.0, Listeria One step CS for surface - 20 tests
	Reveal 2.0, Listeria Enrichment 48 h Complete System

CS = Complete System

1-888-8BIOVET

(824-6838) option 3

order@biovet-inc.com

4375 Beaudry, Saint-Hyacinthe QC J2S 8W2

T. 450 771-7291 • F. 450 771-4158

© Biovet Inc. – June 2022